

SLIT LAMP AMT SL - 01



SPECIFICATION

Microscope	: Galilean
Magnifications Change	: 2 position Lever
Eye pieces	: 10x and 12.5
Magnifications Range	: 10x,16x or 15x,24x
Real field of view in mm	: 18,11 or 15
PD Range	: 54 -90 mm
Diopter Adjustment	: +/-5
Working Distance	: 100 mm
Slit Width	: 0-14mm
Slit Length	: 0-14mm
Slit Apertures	: 0.2, 1,3,4,6,8, 10,14 mm
Slit Angles	: 0° -180°
Slit Inclination	: 5° -10°-15°-20°
Slit Illumination	: 6V-20W
Filters	: Heat absorbing, Green and Cobalt Blue

MOVEMENT RANGES

Longitudinal(In/Out)	: 120mm
Lateral(Left/Right)	: 98 mm
Vertical(Up/Down)	: 56 mm
Voltage	: 110 or 240
Frequency	: 50/60 Hz
Table Base Dimensions	: 522 mm x 342 mm
Weight(packed)	: 26 Kg

KERATOMETER AMT-KM-02

FEATURES

- Coincidence focusing system. Positive measurements to 0.02mm accuracy.
- 15 X Wide-angle piece.
- One position instrument. Measures both meridians without changing optical system.
- Simple vertical adjustable head to fit patient.
- Two-way adjustable head and chin rests.
- Dual eye-level sighting system. Facilitates horizontal alignment.
- Positive fixation. Permits rapid measurement of central corneal area.
- Precision objectives, Assures durability, Achromatic Lenses.
- Measuring range: Radius of curvature, 9.4 ~ 6.4 mm (in 0.05mm steps) Diopter 36~52 (in 0.25 steps).

DIMENSIONS

- H,L,W : 16" x 20" x 7", W : 11 kg.
- Color / Finish : Black or Off White
- Std. Accs. : Dust Cover, 2 - Spare Bulbs



Visit our online store

ophthokart.com
the mega online **Ophthalmic** store



AMT
ALLIED MEDICAL TECHNOLOGIES

Allied Medical Technologies

Plot No. 153 F.I.E Patparganj, Industrial Area, Delhi - 92
Tel./Fax : 0091 - 11 - 42664008
Cell : 0091 - 9811008321, 8586924008, 01, 02
Website : www.amtindia.co, www.ophthokart.com
E-mail : pkrajput2002@yahoo.co.in /
pkrajput2002@gmail.com