



Sport Frames are available in the following finishes



Blue Red Black Sky Blue Pink



Mini Loupes



Standard Loupes
SuperVu Galilean system



Hi Res Loupes
SuperVu Hi-Res Galilean system



Prismatic Loupes
XL-Advantage system

Recommended working distances

Magnification	Working distance	Field of view
2.5x	34cm	7.5cm
	42cm	9.5cm
	46cm	10.5cm
	50cm	11.4cm

Magnification	Working distance	Field of view
2.0x	34cm	11.4cm
	46cm	15.2cm
2.5x	34cm	8.1cm
	42cm	10cm
	46cm	10.9cm
	50cm	11.8cm
3.0x	34cm	5.9cm
	42cm	7.1cm
	46cm	7.7cm
	50cm	8.3cm

Magnification	Working distance	Field of view
3.0x	34cm	6.8cm
	42cm	8.2cm
	46cm	8.8cm
	50cm	9.5cm

Magnification	Working distance	Field of view
3.5x	34cm	5.5cm
	42cm	7.1cm
	46cm	7.9cm
4.5x	34cm	4.6cm
	42cm	6.1cm
	46cm	6.9cm
5.5x	34cm	3.5cm
	42cm	4.7cm
	46cm	5.3cm

Visit our online store

ophthokart.com
the mega online Ophthalmic store



AMT
ALLIED MEDICAL TECHNOLOGIES

Allied Medical Technologies

Plot No. 153 F.I.E Patparganj, Industrial Area, Delhi - 92
Tel./Fax : 0091 - 11 - 42664008
Cell : 0091 - 9811008321, 8586924008, 01, 02
Website : www.amtindia.co, www.ophthokart.com
E-mail : pkrajput2002@yahoo.co.in /
 pkrajput2002@gmail.com