

3 Step / 5 Step Slitlamp



3 Step / 5 Step Slitlamp

- Microscope Type Parallel Galilean
- Magnification changer Revolving the drum
- Eyepiece 12.5X
- Total magnification 10X 16X, 25X / 6x, 10x, 16x, 25x, 40x
- Diameter of Vision field $\phi 23$ $\phi 14$ $\phi 8.7$ (mm) / $\phi 37$, $\phi 23$, $\phi 14$, $\phi 8.7$, $\phi 5.7$
- Pupil adjusting range 54 mm~82 mm
- Diopter adjustment -7D ~ +7D
- Slit image length 0.1mm~14mm
continuous
- Diameter of aperture $\phi 14$ mm, $\phi 10$ mm, $\phi 5$ mm, $\phi 3$ mm, $\phi 1$ mm, $\phi 0.2$ mm
- Slit angle 0~180 Rota table
continuous
- Illumination tilting 5, 10, 15, 20 totally
4 step
- Filter Heat absorption, Grey, red free(Green), Cobalt blue
- Illumination 12v, 30w halogen bulb
- Fixation Red LED 3.5V
- Power Input voltage 100~240V
- Dimension 740*650*450mm
total 25kg

Distributor :



Allied Medical Technologies

Plot No. 153 F.I.E Patparganj, Industrial Area,
Delhi - 110092,
Tel./Fax : 0091 - 11 - 42664008
Cell : 0091 - 9811008321, 8586924008, 01, 02
Website : www.amtindia.co
E-mail : pkrajput2002@yahoo.co.in /
pkrajput2002@gmail.com